

HARTMANN



HydroTac[®] transparent
Transparent hydrogel dressing


SIMPLE AND EFFECTIVE

NEW

HydroTac[®] transparent

Previously named Hydrosorb[®], the dressing's name has now been changed to HydroTac[®] transparent

HydroTac[®] transparent is an individually sealed, sterile wound dressing. Release moisture and absorb a limited amount of exudate. The dressing promotes wound healing by supporting the epithelialization and granulation process.^[1]



Support:

Provides the optimal moisture environment especially for low to non exuding wounds to support wound healing

Soothing:

Provides a soothing and cooling effect to the wound^[3]

Protective:

Dressing's outer layer acts as an effective barrier to bacterial penetration



Skin friendly:

Self-adhesive, with good skin tolerance

Transparent:

Transparent and moisturising dressing

When to use

- For dry or lightly exuding wounds during the granulation and epithelialisation stages
- Chronic and acute wounds – e.g. burns, abrasions, skin tears, etc.

Application

- Easy to handle
- Atraumatic dressing change^[1,2]
- Can remain on the wound up to seven days

HydroTac[®] transparent – available in 4 versions

the non-border version (additional dressing material required for fixation)



5 x 7.5 cm



10 x 10 cm



10 x 20 cm



20 x 20 cm





Helps. Cares. Protects.

Product details

Product	Product size	Product code	PIP code	Pack size
HydroTac® transparent	5 x 7.5 cm	685905	4196481	10
	10 x 10 cm	685906	4196515	10
	10 x 20 cm	685907	4196507	10
	20 x 20 cm	685908	4196499	10



[1] P. Zollner et al (2007) Clinical performance of a Hydrogel dressing in chronic wounds, a prospective observational study. Journal of Wound Care 16, No. 3, March 2007 [2] (2018) Clinical Evaluation Report, Hydrosorb® transparent and hydrosorb® transparent comfort. 1.0/03.07.2018 [3] Bazire et al (2015) Hydrosorb® versus control (water based spray) in the manner of radio-induced skin toxicity: Results of a multi-centre controlled random trial. Radiotherapy and Oncology 117, 229-23

For more information contact your local HARTMANN representative or contact our customer services team

HARTMANN UK & Ireland +44 (0)1706 363200 | Email: info@uk.hartmann.info | www.hartmann.co.uk

Never miss an update from HARTMANN UKI - Join www.hartmannmarketing.com online now

 HARTMANN UKI Healthcare Professionals  HARTMANN GROUP

©2021 Paul Hartmann Limited. ® represents a trademark of Paul Hartmann Limited, a company registered in England and Wales (company number 01523121) WM-000154/21-GB