



HARTMANN



Cosmopor[®] I.V. transparent

Cosmopor[®] I.V. transparent is a non-active, non-invasive and sterile self-adhesive dressing for covering and fixation of cannulae/catheters

Cosmopor® I.V transparent

Skin friendly and easy to apply



New:
fixation strips
included




Visual control
of puncture site
without dressing
change

Remains firmly
in place

Breathable

Barrier against
water and bacterial
penetration



Cosmopor® I.V. transparent is a non-active, non-invasive and sterile self-adhesive dressing for covering and fixation of cannulae/catheters.

When to use

- Covering and fixation of peripheral intravenous cannulae

Product benefits

- High degree of adhesiveness and the included fixation strips allow the dressing to remain firmly in place during the wear time
- Allows visual control of the puncture site without dressing change
- The waterproof characteristics of the film provide a barrier against water, bacterial penetration and external contaminants
- Dressing is breathable and permeable to vapor
- Easy to apply and skin friendly



HARTMANN

Helps. Cares. Protects.


Cosmopor® I.V. transparent



6 x 5 cm



9 x 7 cm



12 x 10 cm

Size**Pack****Hartmann
Code****NHS Code**

6 x 5cm (2.4 x 2")

100 pcs

900817

ELW1019

9 x 7cm (3.5 x 2.8")

100 pcs

900818

ELW1020

12 x 10cm (4.7 x 4")





50 pcs

900819

ELW1022






Cosmopor I.V.®

6 x 5cm | 9 x 7cm



Cosmopor I.V.®

12 x 10cm



For more information contact your local HARTMANN representative or contact our customer services team

HARTMANN UK & Ireland +44 (0)1706 363200 | Email: info@uk.hartmann.info | www.hartmann.co.uk

Never miss an update from HARTMANN UKI - Join www.hartmannmarketing.com online now

HARTMANN UKI Healthcare Professionals HARTMANN GROUP

©2022 Paul Hartmann Limited. ® represents a trademark of Paul Hartmann Limited, a company registered in England and Wales (company number 01523121) WM-000071/22-GB