

HARTMANN



HydroTac[®] transparent
Transparent hydrogel dressing


SIMPLE AND EFFECTIVE

NEW

HydroTac[®] transparent

HydroTac[®] transparent comfort

HydroTac[®] transparent and transparent comfort are individually sealed, sterile wound dressings. Both versions release moisture and absorb a limited amount of exudate. The dressings promote wound healing by supporting the epithelialization and granulation process.^[1]



Support:

Provides the optimal moisture environment especially for low to non exuding wounds to support wound healing

Soothing:

Provides a soothing and cooling effect to the wound^[3]

Protective:

Dressing's outer layer acts as an effective barrier to bacterial penetration



Skin friendly:

Self-adhesive, with good skin tolerance

Transparent:

Transparent and moisturising dressing

When to use

- For dry or lightly exuding wounds during the granulation and epithelialisation stages
- Chronic and acute wounds – e.g. burns, abrasions, skin tears, etc.

Application

- Easy to handle
- Atraumatic dressing change^[1,2]
- Can remain on the wound up to seven days

HydroTac® transparent – available in 4 versions the non-border version (additional dressing material required for fixation)



5 x 7.5 cm



10 x 10 cm



10 x 20 cm



20 x 20 cm

HydroTac® transparent comfort – available in 4 versions features an adhesive border, coated with an acrylic based adhesive. Thus, no secondary dressing is required for fixation.



8 x 8 cm



12.5 x 12.5 cm



10 x 20 cm



6.5 x 10 cm





Helps. Cares. Protects.

Product details

Product	Product Sizes		Product code	NHSSC code	PIP code	Pack size
	Overall	Wound pad size				
HydroTac [®] transparent	5 x 7.5 cm	n/a	6859050	ELE109	4196481	10
	10 x 10 cm	n/a	6859060	ELE110	4196515	10
	10 x 20 cm	n/a	6859070	ELE106	4196507	10
	20 x 20 cm	n/a	6859080	ELE111	4196499	10
HydroTac [®] transparent comfort	8 x 8 cm	4 x 4 cm	6859250	ELE112	-	10
	12.5 x 12.5 cm	8 x 8 cm	6859260	ELE107	-	10
	10 x 20 cm	5.5 x 15.5 cm	6859270	ELE108	-	10
	6.5 x 10 cm	2.5 x 6 cm	6859280	ELE118	-	10



[1] P. Zollner et al (2007) Clinical performance of a Hydrogel dressing in chronic wounds, a prospective observational study. Journal of Wound Care 16, No. 3, March 2007 [2] (2018) Clinical Evaluation Report, Hydrosorb[®] transparent and hydrosorb[®] transparent comfort. 1.0/03.07.2018 [3] Bazire et al (2015) Hydrosorb[®] versus control (water based spray) in the manner of radio-induced skin toxicity: Results of a multi-centre controlled random trial. Radiotherapy and Oncology 117, 229-23

For more information contact your local HARTMANN representative or contact our customer services team

HARTMANN UK & Ireland +44 (0)1706 363200 | Email: info@uk.hartmann.info | www.hartmann.co.uk

Never miss an update from HARTMANN UKI - Join www.hartmannmarketing.com online now

HARTMANN UKI Healthcare Professionals HARTMANN GROUP

©2022 Paul Hartmann Limited. ® represents a trademark of Paul Hartmann Limited, a company registered in England and Wales (company number 01523121) WM-000111/21-GB