

HARTMANN



Hydrofilm®

Your discreet and waterproof transparent film



Hydrofilm®

Sterile film barrier to reduce the risk of infection and protect against physical damage on dry and primary healing wounds.


- **Transparent** - Allows assessment of surrounding skin without removal of dressing
- **Impermeable to microorganisms** - Viruses, bacteria and external contaminants
- **High breathability** - Facilitates the free passage of vapours and gases



Provides protection during shower

Almost painless removal by simply overstretching the film

4-step application



Size	Packing contents	REF	PIP codes	NHSC codes
6 x 7 cm	10	970000	3426665	
6 x 7 cm	100	970001	4326665	
10 x 12.5 cm	10	970002	3426228	
10 x 12.5 cm	100	970003		
10 x 15 cm	10	970004	2667350	
10 x 15 cm	50	970005		
15 x 20 cm	10	970006	3426244	
15 x 20 cm	50	970007	3426244	
10 x 25 cm	25	970008	3426236	
12 x 25 cm	25	970009	3427368	
20 x 30 cm	10	970010	3426251	ELW248

For more information contact your local HARTMANN representative or contact our customer services team

HARTMANN UK & Ireland +44 (0)1706 363200 | Email: info@uk.hartmann.info | www.hartmann.co.uk
 Never miss an update from HARTMANN UKI - Join www.hartmannmarketing.com online now

HARTMANN UKI Healthcare Professionals HARTMANN GROUP

©2023 Paul Hartmann Limited. ® represents a trademark of Paul Hartmann Limited, a company registered in England and Wales (company number 01523121) WM-000018/23-GB March 2023

